



Transfusion

Analytical Performance Specifications



Foetal Maternal Haemoglobin (FMH)

Measurand	Lower Limit	Upper Limit
Kleihauer %	$\pm 0.06 \leq 0.18 \%$	$\pm 40\%$ if $> 0.18\%$
Kleihauer mL	$\pm 1.4 \leq 4.3 \text{ mL}$	$\pm 40\%$ if $> 4.3 \text{ mL}$
Kleihauer VIAL	$\pm 1.0 \leq 9.0$	$\pm 10\%$ if > 9.0
Flow %	$\pm 0.05 \leq 0.20 \%$	$\pm 25\%$ if $> 0.20 \%$
Flow mL	$\pm 1.1 \leq 4.4 \text{ mL}$	$\pm 25\%$ if $> 4.4 \text{ mL}$
Flow VIAL	$\pm 1.0 \leq 9.0$	$\pm 10\%$ if > 9.0

